I

Italiano

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TERMOMETRO

DIGITALE

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DIGITAL

THERMOMETER

Before use, read the instructions and warnings given in this handbook carefully and keep the handbook in a safe place.

Digital thermometer with a large LCD display for easy reading. It measures body temperature in 10 seconds. Fast and accurate fitted with a sound signal when the temperature reading has been taken. Stores the last read temperature, to evaluate any variations. Extremely easy to use, it is kept in a practical case to avoid risk of breaking. The probes can be washed, and it turns of automatically.

■ GENERAL WARNINGS

- 1 Self-measuring means control, not diagnosis and treatment, any unusual readings must be discussed with your doctor. Under no circumstances must you change the medicine doses that have been prescribed by your doctor.
- 2 Do not squeeze, bend, drop or break the thermometer.
- 3 Keep it out of the reach of children.
- 4 Do not expose to sunrays or direct sources of heat. Keep in a dry clean place.
- 5 Do not immerse the thermometer in water or other liquids, only the probe is waterproof.
- 6 This thermometer reads only the armpit, rectal or oral temperature. Do not use it for reading auricular temperatures.
- 7 Do not drink hot or cold drinks, do not smoke, do not take a shower before measuring body temperature.

■ INSTRUCTIONS FOR USE
* IN THE MOUTH AND IN THE RECTAL

- 1 Press the ON/OFF key. For a few seconds all the function symbols light up on the display. After the thermometer displays 36.5°C and L. When the °C symbol starts flashing, the thermometer is ready for use.
- 2 Place the probe as follows:
IN THE MOUTH: place the probe under the tongue and keep the mouth closed for about 10 seconds, until the temperature stabilises. The standard temperature for this type of reading is 37°C.
UNDER IN THE RECTAL: lubricate the probe with a water-soluble gelatinous solution. Do not use Vaseline. Insert the probe into the rectal for about 1 cm., until the temperature stabilises (about 10 seconds). The standard temperature for this type of reading is 37.0°C.

* UNDER THE ARM

To measure temperature in the armpit area, proceed as follows:

- 1 Dry the armpit well. Place the probe under the armpit and keep the arm well tightened to the chest. (In this way, room temperature does not influence the correct measurement of body temperature).
- 2 Press the ON/OFF key. All function symbols light on the display, just for few seconds. Then, the thermometer displays 36.5°C and L. When the symbol °C blinks, the device is ready for use. Temperature stabilizes after about 10 seconds. Standard temperature for this measurement type is 36.7°C.

- 3 When the thermometer gives off the sound signal, remove it. The reading is terminated and shown on the display. The thermometer turns off automatically after 10 minutes without use (in the meantime the reading is still displayed) otherwise you can turn it off by pressing the ON/OFF button. If the temperature were above 44°C, the display shows HC; if it were below 32°C, the temperature shows LC.

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■ MEMORY FUNCTION

The thermometer stores the last reading in its memory, to evaluate any variations. To recall the last reading, press the ON/OFF button for a few seconds until the reading appears on the display; in the bottom right corner of the display, "M" appears. Release the ON/OFF button to continue with measuring the actual temperature.

■ MAINTENANCE

Clean the thermometer probe with a cloth dampened with a 0.02% hypochlorite sodium solution or with alcohol or with any type of topical disinfectant before and/or after each use. The rest of the thermometer must be cleaned with a dry cloth. Keep the thermometer safe in its special case.

■ TECHNICAL DATA

Measuring range: 32°C-44°C
Tolerance: +/- 0.1°C 35°C-39.0°C, otherwise +/- 0.2°C
Accuracy: 0.1°C
Working conditions: 0°C to 50°C, Humidity: 10-83%
Conservation conditions: -20°C to 60°C, Humidity: 10-83%
Battery autonomy: approx. 100 hours continuous use
Battery consumption: 0.15 mW
Display: liquid crystal, 3 1/2 characters.
Size: 126.5 mm x 18 mm x 10 mm
Weight: 7 g (with battery)

■DISPOSAL PROCEDURE

 The symbol located on the bottom of the device indicates the separated collection of electric and electronic equipment. At the end of the useful life of the device, do not dispose it as mixed solid municipal waste, but dispose it in a specific collection centre located in your area. This procedure of separated collection of electric and electronic equipment is carried out foreseeing a common EU environmental policy aiming at safeguarding, protecting and improving the quality of environment.

■WARRANTY

This appliance is warranted for a period of 2 years from the date of purchase as shown by the **dealer's stamp and signature and by the purchaser's receipt for tax purposes which must be kept attached to this warranty**. The warranty covers only any faults in manufacturing and does not cover accidental damage, incorrect use, negligence or misuse of the product. Use only the accessories provided; the use of accessories other than the approved ones can make your warranty void. Never open the appliance; in the event of any opening or tampering the warranty will automatically become invalid. The warranty does not cover either the parts subject to wear and tear or the batteries (if provided). The warranty expires 2 years from the date of purchase; after this period any technical service shall be charged. Information regarding technical service either under the terms of this warranty or to be paid for can be obtained by contacting info@laica.com. No compensation of any kind is owed for repairs or substitutions of products under the terms of this guarantee. In the event of faults, go back to the retailer; DO NOT send the appliance directly to LAICA. Use in public places is not allowed. The manufacturer shall not be liable for any damages to property, property or pets, either directly or indirectly, caused by the non compliance with the instructions laid down in the instruction manual, in particular the warnings regarding the installation, use and maintenance of the appliance. It is Laica's policy to improve its products constantly; Laica therefore reserves the right to change its products, either partially or totally, without previous notification and according to the manufacturing requirements, without either Laica or its retailers being in any way liable for it.

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è un marchio LAICA SpA